

MYCOTOXINS

AgraStrip® Pro WATEX® AgraVision™ Pro



The AgraStrip® Pro WATEX® test system enables the rapid and simple on-site quantification of mycotoxins in a variety of agricultural commodities. Both extraction and assay have been streamlined, reducing steps to a bare minimum while setting new standards in accuracy and usability.



Product Description

AgraStrip® Pro WATEX® test kits are lateral flow devices for the rapid, on-site detection of mycotoxins in a variety of agricultural commodities.

Strips are available for total aflatoxins (B1, B2, G1, G2), deoxynivalenol, zearalenone and total fumonisins (B1, B2, B3).

With the newly developed AgraVision™ Pro reader, the system features a shared extraction and dilution procedure and an intuitive walk-away operation to allow simple testing for several different mycotoxins in less than 10 minutes while helping to eliminate common handling errors.

The system can independently and simultaneously test up to 4 samples. The assay takes only 4 minutes, and the reader controls timing, temperature and flow, all in a streamlined procedure. With its integrated incubator and resistance to dust and dirt, the reader secures precise and quantitative results and allows consistent documentation.

Our proprietary water-based extraction eliminates hazards both to the user and the environment.

Applications

A variety of commodities can be analyzed. Please refer to current application briefs for your commodity and mycotoxin of interest.

Kit Contents

AgraStrip® Pro WATEX® test kits come complete with strips, strip cartridges, microcentrifuge tubes and a holder, a dilution buffer, extraction buffer bags, pipette tips and Whirl-Pak® filter bags to analyze up to 40 samples.

FEATURES & BENEFITS

Rapid

Up your grain intake with quick sample preparation, a 4-minute assay time and the capacity to test 4 samples at once on-site.

Simple

Eliminate errors with a common extraction/dilution and an intuitive, walk-away operation. No sophisticated training process required. The water-based extraction also means no handling or disposal of organic solvents.

Robust

Get precise, reliable results in even the toughest of conditions thanks to the integrated incubator and resistance to dust and dirt.

Connected

Stay on top of your data with blockchain-enabled test strips. Transfer your data to your PC or print it with ease.

Outstanding Support

You're not alone. Whether it's for validations, questions or custom solutions, our experts are there, whenever you need them.

AgraStrip® Pro WATEX® AgraVision™ Pro



WATEX®
One for All
WATER EXTRACTION

To ensure accurate results, please follow the package insert that came with your test kit.

3-Step Procedure in Brief

1 Extract your sample



Add sample and buffer bag to the Whirl-Pak® filter.



Add distilled, deionized or bottled water.




Shake and let settle.

4
minute
assay

2 Perform the test-specific dilution and centrifugation



Add assay dilution buffer to a microcentrifuge tube.



Transfer sample extract into the microcentrifuge tube and centrifuge for 30 seconds.

3 Run the assay with the AgraVision™ Pro reader



Transfer the cartridge with the strip to the reader and enter sample and test information.



Add sample extract to the cartridge. The reader starts the assay and displays results when finished.

Ordering Information - AgraStrip® Pro WATEX®/AgraVision™ Pro

Product	Description	Item Number	Limit of Detection	Quantitation Range
AgraStrip® Pro Total Aflatoxin WATEX®	40 tests	10006322	2/1 ppb*	0 – 460 ppb
AgraStrip® Pro Deoxynivalenol WATEX®	40 tests	10006320	0.1 ppm	0 – 44 ppm
AgraStrip® Pro Total Fumonisin WATEX®	40 tests	10006318	0.1 ppm	0 – 44 ppm
AgraStrip® Pro Zearalenone WATEX®	40 tests	10006316	25 ppb	0 – 1650 ppb
AgraVision™ Pro Reader	Reader	10006324		n/a
AgraVision™ Pro Printer	Printer	10006483		n/a
AgraVision™ Pro Scanner	Scanner	10006484		n/a

* High-sensitive method for aflatoxin available

Romer Labs Inc · Newark, DE, USA · T +1 302 781 6400 · F +1 302 781 6378 · E office@romerlabs.com

Romer Labs Diagnostic GmbH · Tulln, Austria · T +43 2272 615 33 · F +43 2272 615 33 13177 · E office-europe@romerlabs.com

Romer Labs Singapore Pte. Ltd · Singapore · T +65 6631 8018 · F +65 6275 5584 · E salesasia@romerlabs.com

www.romerlabs.com